

**Product Name :**  
Submersible Pump Test Rig

**Product Code :**  
TOURBO SUK-AS03LAB003



**Description :**

Submersible Pump Test Rig

**Technical Specification :**

Submersible Pump Test Rig: The designed and developed submersible pump Test Rig is a self-contained unit operated on closed circuit basis containing a sump tank. The set-up consists of a submersible pump coupled with a DC Motor. Power input to the DC Motor is varied by means of a Thruster controlled DC Drive to vary the RPM of motor. A RPM Indicator with Proximity sensor indicates the RPM of Pump. Flow of water is measured by using measuring tank and stopwatch. Vacuum gauge



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265