

Product Name :
Computerized Axial Flow Pump

Product Code :
TOURBO SUK-AS03LAB005



Description :

Computerized Axial Flow Pump

Technical Specification :

Computerized Axial Flow Pump: Fluid Machinery Lab equipment: Turbo Machinery lab equipment: In an axial-flow pump the pumping medium flows through the impeller (here a propeller) in the axial direction. In axial-flow pumps, the pressure build-up is not achieved by the centrifugal force but like the aerodynamic principle, by the propeller blade. Therefore these pumps are also known as propeller pumps.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265