

Product Name :
Computerized Kaplan Turbine Test Rig

Product Code :
TOURBO SUK-AS03LAB008



Description :

Computerized Kaplan Turbine Test Rig

Technical Specification :

Computerized Kaplan Turbine Test Rig: Fluid machinery lab equipment :
Tourbo machinery lab equipment : The designed and developed experimental setup helps to investigate the characteristic behavior of a simple-regulated Kaplan turbine during operation. The trainer is provided with a closed water circuit with tank, submersible pump and throttle valve for adjusting the flow rate



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265