

**Product Name :**  
Computerized Reciprocating Pump Test Apparatus

**Product Code :**  
TOURBO SUK-AS03LAB009



**Description :**

Computerized Reciprocating Pump Test Apparatus

**Technical Specification :**

Computerized Reciprocating Pump Test Apparatus: Fluid machinery lab equipment : Turbo machinery lab equipment : The test rig consists of a closed water circuit with water tank, a piston pumps with variable speed via a frequency converter and an air vessel. The piston of the pump is mounted in a transparent housing and can be observed during operation. The cycle that takes



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265