

Product Name : UV Light Testing Meter	Product Code : MPI SUK-AS03LAB005
---	---



Description :

UV Light Testing Meter

Technical Specification :

Application

Weight (with sensor)

Dimensions 25 mm

Battery giving 60 - 100 Operating Hours

Operating Temperature

Central Wavelength

Batteria voltaic Cell with UV pass filter

Accuracy

Display

352Mm* 66mm* 33mm

Usage Light Meter, Lux Meter



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265