

Product Name :
Digital Temperature Scanner With Logger

Product Code :
PROCESS SUK-AS03LAB003



Description :

Digital Temperature Scanner With Logger

Technical Specification :

~~Temperature Range~~

~~Weight~~ 1500 Gms

~~Size~~ 96 X 120 Mm

~~Electric Source~~

~~Digital Type~~

~~Power~~ - 230 V



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265