

**Product Name :**  
Hot Specimen Mounting Press

**Product Code :**  
METALLURGICAL SUK-  
AS03LAB007



**Description :**

Hot Specimen Mounting Press

**Technical Specification :**

Dimensions 150x150x500 Cm

Voltage 220V

Temperature Setting Range 300-1000°C

Net Weight 16 Kg

Max Power 1800W



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265