

Product Name :
Digital Count Tester

Product Code :
TESTING SUK-AS03LAB001



Description :
Digital Count Tester

Technical Specification :

Overview
As one of the most promising manufacturer and supplier in the current market scenario, our organization is capable of providing a wide range of Digital Count Tester. This is basically designed for counting the unit counts in any basic and measure the count accurately without any error. This device is used in numerous applications. Our offered Digital Count Tester is manufactured in accordance with industry quality standards by making use of high grade components along with modern machine and high end techniques.

Techniques:

Features:

- Excellent working
- Dimensional stability

• Light weight



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265