

**Product Name :**  
Digital Count Tester

**Product Code :**  
TESTING SUK-AS03LAB001



**Description :**

Digital Count Tester

**Technical Specification :**

Overview  
As one of the most promising manufacturer and supplier in the current market, our organization is capable of providing a wide range of Digital Count Tester. This is basically designed for counting the unit counts in any form and measure the count accurately. The device offers excellent accuracy and precision. The device offers excellent accuracy per count. This device used in numerous fields and industries. Our latest Digital Count Tester is manufactured in accordance with industry quality standards by making use of high grade components along with modern machine and high end techniques.

Techniques:

Features:

- Excellent working
- Dimensional stability

• Light weight



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265