

Product Name :
Manual Crockmeter

Product Code :
TESTING SUK-AS03LAB002



Description :

Manual Crockmeter

Technical Specification :

Manual Crockmeter

Crock
Front Operated, provides a reciprocating rubbing motion simulating the action of a human finger and thumb. It is used for determination of loss of fluid and dried matter in dry and wet rubbing. Crock
Front Operated is designed for the 10 cm diameter top cover back and forth, with each complete turn of the crank, in a straight line along a 100 mm x 50 mm on the specification, with a downward force of

2N




Accessories

- Main Unit - 01 No.
- Cooking Cloth (200 PCS PACK)
- S.S Spring Clips (For mounting Crock square) - 02 Nos.
- Blotting papers (200 PCS PACK)
- Specimen Holder - 01 No.

- User's Manual : 01 No.

- Calibration Certificate : 01 No.

Technical Specifications



Color (if available)

Site or Test Specimen

Depth/Altitude of the Site

Depth Weight of the Site

Depth/Altitude of the Site

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com