

Product Name :
Digital Crockmeter

Product Code :
TESTING SUK-AS03LAB003



Description :

Digital Crockmeter

Technical Specification :

Giving precedence to the varied requirements of the industry, we have been successfully engaged in manufacturing and supplying an extensive collection of Digital Crockmeter. This has been designed for offering a reciprocating rubbing motion that simulates the action of finger and thumb. Therefore, this product is used for determining hardness of glass and pressed bodies against dry or wet rubbing. Our crockmeter is designed under the close guidance of skilled professionals by making use of excellent grade raw material.

Our clients can use this Digital Crockmeter in numerous specifications at varied prices.


Feature:

- Easy usability
- Light weight
- Simple to read




Accessories

- Main Unit: 01 No.
- Counting Cheete (200 PCS PACK)
- S.S Spring Clips (For mounting Clock square): 02 Nos.
- Blotting papers (200 PCS PACK)
- Specimen Holder : 01 No.
- User's Manual : 01 No.
- Calibration Certificate (Traceable to NPL) : 01 No.

Technical specification



Diameter of the Rubbing Finger



Packing Dimensions of the box

Drop Weight of the box

Drop/Minic/Maxic Box



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265