

Product Name :
Full Shoe Flex Tester

Product Code :
FOOTWEAR SUK-AS03LAB001



Description :

Full Shoe Flex Tester

Technical Specification :

Being one of the distinguished names in the industry, we are engaged in manufacturing, supplying and exporting Full Shoe Flex Tester.

Our range is manufactured using superior quality of raw materials and components. These equipments are broadly used to measure the ability of the full shoe to withstand the effect of flexing stresses produced on the various parts of the shoe. Acknowledged for their durability, high accuracy and long service life, our offered Full Shoe Flex Tester is high in demand across the market.

Features:

Although the various parts like uppers and soles of a shoe are tested separately for flexing endurance with the help of flexometer and Ross flex tester, the shoe flex tester gives results which are more related to the actual

performance of the shoe as a whole.

Technical specification

60 Digit
Electric Source
50 Material
230 Voltage
50 Frequency
Blade Enterprises
80 Maximum Eccentricity

55 Deg (Max)
20 Cycles/Minute
Digital Type



Astra Scientific

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265