

**Product Name :**  
Full Shoe Flex Tester

**Product Code :**  
FOOTWEAR SUK-AS03LAB001



**Description :**

Full Shoe Flex Tester

**Technical Specification :**

Being one of the distinguished names in the industry, we are engaged in manufacturing, supplying and exporting Full Shoe Flex Tester.

Our range is manufactured using superior quality of raw materials and components. These equipments are broadly used to measure the ability of the full shoe to withstand the effect of flexing stresses produced on the various parts of the shoe. Acknowledged for their durability, high accuracy and long service life, our offered Full Shoe Flex Tester is high in demand across the market.

**Features:**

Although the various parts like uppers and soles of a shoe are tested separately for flexing endurance with the help of flexometer and Ross flex tester, the shoe flex tester gives results which are more related to the actual

performance of the shoe as a whole.

Technical specification

600g
Electric Source
Material
230V
50Hz
Blade Enterprises
Maximum Eccentricity

55 Deg/Sec (Max Angle)
200 Cycles/Minute
Digital Type



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265