

| | |
|---|--|
| Product Name : Ross Flex Tester | Product Code : FOOTWEAR SUK-AS03LAB004 |
|---|--|



| |
|----------------------|
| Description : |
| Ross Flex Tester |

| |
|--|
| Technical Specification : |
| The ability of shoes or chappals to withstand repeated bending experience during walking without developing a crack within their estimated life span is a very important factor for determining the quality of its sole. Also, the ability of the sole to withstand propagation of any slit formed due to a sharp edge or pointed object during walking is equally important in assessing its quality. |

| | | |
|--|--------------------|--|
| Technical specification | | |
| <table><tr><td>Bend Radius</td></tr><tr><td> </td></tr></table> | Bend Radius | |
| Bend Radius | | |
| | | |

| |
|------------------------------|
| 220 Age |
| 500 Hz |
| 150 x 25 x 6.25 mm Diameter |
| 3.5 Width |
| Flexing 5 Cycles Per Minutes |
| 90 Degree Specimen Angle |



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265