

Product Name :
Sole Adhesion Tester

Product Code :
FOOTWEAR SUK-AS03LAB005



Description :

Sole Adhesion Tester

Technical Specification :

Sole Adhesion Tester.

Under the strict supervision of our well trained work force, these instruments are manufactured using supreme quality materials. These are widely used in shoe factories or at lab to determine the adhesion strength of sole at the toe and heel. **Sole Adhesion Tester** is a time saving tools and highly demanded for their accuracy, long functional life and outstanding performance.

Further details:

The Strength Of Adhesion Of A Sole To Its Upper Is Determined By Measuring The Force Needed To Pull Off The Sole From The Upper. A Suitable Last Is Placed Inside The Shoe, Which Is Then Placed On A Straight Edge Acting As A Fulcrum For Applying Load On The Joint Of The Sole And Upper. The Maximum Load The Sole Can Withstand Before Coming Off Is Determined And Reported As The Sole Adhesion Strength.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265