

<b>Product Name :</b> Bennewart Flex Tester	<b>Product Code :</b> FOOTWEAR SUK-AS03LAB007
--	--



<b>Description :</b>
Bennewart Flex Tester
<b>Technical Specification :</b>
<p><b>Bennewart Flex Tester.</b></p> <p>These equipments are manufactured by our highly skilled manpower using sophisticated technology. Our quality analyzer stringently checks these flex testing equipments on quality parameters before delivering to our clients. <b>Bennewart Flex Tester</b> is broadly used to measure the resistance of a component to cut growth during repeated flexing. Clients can buy these equipments from us t cost effective prices.</p> <p><b>Further details:</b></p> <p>It can also be used to assess the effect of surface patterns on crack initiation and growth. This test is especially applicable to the outsoles of footwear, but may also be used with certain other flexible components.</p>

Technical specification

<b>Bending Mandrel Diameter</b>
<b>Display Type</b>
<b>Material</b>
<b>Grade</b> Automatic
<b>Voltage</b>
<b>Frequency</b>
<b>Brand</b> Enterprises
<b>Operation Frequency</b>

Flare Angle	
Digitizer	



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265