

<b>Product Name :</b> Flexometer Machine	<b>Product Code :</b> FOOTWEAR SUK-AS03LAB008
---	--



**Description :**

Flexometer Machine

**Technical Specification :**

In this test, test specimens in shape of rectangular pieces are folded and clamped at each end to maintain them in a folded position in a set of grips, one of which is fixed while the other is able to oscillate. The movement of the oscillating grip causes the fold in the test specimen to run along its centre. This operation is carried out repeatedly and the test specimens inspected periodically to assess the damage produced. The equipment used for this test is called as Flexometer.

**Technical specification**

Digital Type
--------------

230V

Voltage

50Hz

Frequency

Blade Enterprises

Test Specimen Dimension

Blade Phase

0.5 to 12.5 Degree

Angle



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265