

Product Name :
Heat Deflection Tester

Product Code :
PLASTIC SUK-AS03LAB005



Description :

Heat Deflection Tester

Technical Specification :

- IS 2267 1972 Specification for Polystyrene Moulding Materials.
- ASTM D 1238 x 1982 Standard Test Method for Flow Rates of Thermoplastics by Extrusion Plastometer.

Technical specification

Dimensions of the test specimen
Diameter of the test specimen, 120 ± 10 mm long

Distance between supporting edges

Radius of supporting and loading edges

Temp 3000 (least count 1°C)



Astra Scientific

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265