

**Product Name :**  
Heat Deflection Tester

**Product Code :**  
PLASTIC SUK-AS03LAB005



**Description :**

Heat Deflection Tester

**Technical Specification :**

- IS 2267 1972 Specification for Polystyrene Moulding Materials.
- ASTM D 1238 x 1982 Standard Test Method for Flow Rates of Thermoplastics by Extrusion Plastometer.

**Technical specification**

Dimensions of the test specimen Diameter of the test specimen, 120 ± 10 mm long
Distance between supporting edges 100 mm
Radius of supporting and loading edges R = 6.25 mm

Temp 3000 (least count 1°C)



Astra Scientific

## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265