

**Product Name :**  
Bottle Compression Tester

**Product Code :**  
PLASTIC SUK-AS03LAB007



**Description :**

Bottle Compression Tester

**Technical Specification :**

With an aim to be a market leader in our domain, we are growing as one of the eminent companies engaged in manufacturing and exporting Bottle Compression Tester. Used in plastic bottles making industries, it is used to test the amount of pressure required to compress thermoplastic bottle. This tester is used by placing specimen between a compression plate and base. Then pressure is slowly applied and readings are recorded by a digital screen. Bottle Compression Tester determines ability of bottles to bear stacking load and impact during transit.

**Features:**

- Can also be used for cartons or paper cups
- Ergonomic structure
- Long lasting life



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265