

Product Name :
Round Sample Cutter

Product Code :
FABRIC SUK-AS03LAB002



Description :

Round Sample Cutter

Technical Specification :

The Round Sample Cutter cuts out circular test specimens of 100 cm² area rapidly and accurately. Grammage or GSM of the material under test is determined by multiplying the mass of the circular specimen obtained by 100.

The equipment consists of an aluminium body fitted with a set of four replaceable blades and also provided with safety catch lock. The reversible blades allow the use of all the four edges of each blade. Two special cutting pads are provided with the unit.

The sample to be cut is placed on the cutting pad on which the sample cutter is placed. After releasing the safety catch lock, a slight pressure is exerted on the hand wheel to bring the blades into contact with the fabric. The hand wheel is then rotated under pressure to cut out the specimens

Technical details

<div>Ocular Specimen Area</div>
<div>Frequency</div>
<div>Material</div>
<div>Voltage</div>



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265