

Product Name : Martindale Abrasion Tester	Product Code : FABRIC SUK-AS03LAB017
---	--



Description :

Martindale Abrasion Tester

Technical Specification :

- Material
Metal
- Power Source
Electric
- Weight
25-40 Kilograms (kg)
- Usage

resistance to abrasion by calculating the loss in mass of the fabric

We have carved a reputed niche in industry for manufacturing and exporting **Martindale Abrasion Tester**. It is most commonly utilized in textiles and leather processing industries. This tester is required to determine wear resistance of a specific fabric or material. It is ideal to be used for testing upholstery, footwear, & gloves as well. **Martindale Abrasion Tester** uses standardized abrasants against which specimens are rubbed to determine their resistance. It is highly appreciated for allowing easy switching between abrasion, pilling, & crocking test tools.

Features:

- Assures quick motion change without need of any tool
- Equipped with LCD display to show results
- User friendly operation

Martindale Abrasion Tester

Technical Specifications

Maximum number of test specimens that can be tested at a time
Speed of outer pegs
Ratio of speeds of outer pegs to central peg
Stroke of pegs
Working area of specimen holder
Pressure on specimens

Related Standards	
ISO 12673 and ISO 12673-1 Abrasion Resistance - Methods for Determination	
BS 5969-1 Method of Test for Abrasion Resistance of Fabrics	



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265