

Product Name :
Pressure Head Tester

Product Code :
FABRIC SUK-AS03LAB018



Description :

Pressure Head Tester

Technical Specification :

- Specimen Size

200 mm Millimeter (mm)

- Material
Metal
 - Power Source
Electric
 - Weight
20-40 Kilograms (kg)
 - Usage
Head Pressure Testing
- IS 1322: 1982 Specification for Bitumen Felts for Water Proofing and Damp Proofing.

Technical details

Pressure Head
Assumed Water Column
Set on Auto Grade
Usage

50 Hz Frequency
200 mm Specimen Diameter
70 s Duration
100 mm Area Diameter



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265