

Product Name :
Traverse Thread Counter

Product Code :
FOAM SUK-AS03LAB002



Description :

Traverse Thread Counter

Technical Specification :

The linear number of cells in any sample of polyurethane foam is determined by finding the average number of cells between two lines spaced a known distance from each other in three mutually perpendicular planes and finding their mean.

The volumetric number of cells in any sample of polyurethane foam is determined by finding the average number of cells between two lines spaced a known distance from each other in three mutually perpendicular planes and multiplying them with each other.

Technical specification

--

50 Hz
Material
Digital Type
Voltage



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265