

Product Name :
Pounding Tester

Product Code :
FOAM SUK-AS03LAB006



Description :

Pounding Tester

Technical Specification :

- Dimensions of test specimen: 380 x 380 x 50 mm
- Diameter of indenter : 250 ± 1 mm
- Number of specimens: One
- Radius at edge of indenter: 25 ± 1 mm
- Movement of indenter : Adjustable between 20 and 50 mm
- Mass of indenter assembly : 76.45 ± 2 kg (force applied - 750 ± 20 N)
- Pounding frequency : 70 ± 2 cpm
- Motor : $\frac{1}{2}$ HP single-phase 230 volts AC



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265