

**Product Name :**  
Tensile Strength Tester

**Product Code :**  
RUBBER SUK-AS03LAB001



**Description :**

Tensile Strength Tester

**Technical Specification :**

We have achieved high acclamation in the industry for offering a top-of-the-line gamut of Tensile Strength Tester.

Manufactured as per industrial standards using quality components, these testing equipments are known for their high efficiency and accuracy. We are offering these testing equipments to our clients without any delay. Our offered Tensile Strength Tester is easy to use and widely used to determine the elongation & tensile strength of variety of raw materials such as fabrics, wires, leather, plastics, belts and rubbers.

Further details:

Tensile strength and elongation are the two prime characteristics of most of the raw materials, . These two properties often play a major role in determining the suitability of any raw material for any specified

application. It is, therefore, of utmost importance to determine these characteristics accurately, conveniently, and quickly.

Technical specification

Stainless Steel
20 Age
50 Hz



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265