

**Product Name :**  
Tension Set Apparatus

**Product Code :**  
RUBBER SUK-AS03LAB003



**Description :**

Tension Set Apparatus

**Technical Specification :**

- Material

SS

- Weight  
5 Kilograms (kg)
  
- Usage  
Rubber

With a vision to be a pioneer in our domain, we have been engaged in manufacturing and supplying Tension Set Apparatus.

It is used for testing tensile strength of various elastomers including rubber, belts, films, and engineered plastics as well. This apparatus is provided with enlarged ends to assure perfect gripping of specimen. Tension Set Apparatus is comprised of a load cell, piston, cylinder, crosshead, notched column, and tension crosshead.

Features:

- Assures correct results
- Easy to use
- Requires low maintenance

Technical specification

Material
Usage
Power Source

Number of Specimen Tested at One Time
Maximum Separation of Grip
Frequency



Astra Scientific

## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265