

Product Name :
Specimen Cutting Press

Product Code :
RUBBER SUK-AS03LAB004



Description :

Specimen Cutting Press

Technical Specification :

Persistently marching in competitive market, our company is recognized as one of the well established manufacturers, suppliers and exporters of Specimen Cutting Press.

Manufactured using sophisticated technology; this cutting press is known for its robust construction, high operational fluency and easy maintenance. We are committed to offer damage free and flawless quality cutting press at the clients' end. Specimen Cutting Press is extensively used for cutting wire, plastic, rubber, leather, plywood and hardwood.

Further details:

Test

specimens of specified shapes and sizes are often needed for conducting various physical tests on a wide variety of materials. These test specimens are normally cut out of the samples of materials under test using appropriate means.

Test

specimens of rubber, leathers, fabrics, foam, flexible or semi-rigid plastics and similar other materials are usually cut with the help of knife edged dies of appropriate shapes and dimensions. The test sample / slab is placed over a flat and smooth piece of , acrylic, or hard rubber (ebonite), while the cutting die is pressed down from the top to cut out the test specimens.

Technical specification

20 Age
Pressure Gauge For Testing Upto 160 kg
1000 kg Compression Load
Material
50 Hz



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265