

**Product Name :**  
Compression Set Apparatus(Quick Clamping Type)

**Product Code :**  
RUBBER SUK-AS03LAB005



**Description :**

Compression Set Apparatus(Quick Clamping Type)

**Technical Specification :**

It required to determine the ability of elastomeric products to maintain elastic properties when kept in compressive state for a long period of time. This apparatus is used by keeping specimen between compression plates and applying desired amount of pressure by turning provided knob. Compression Set Apparatus is ideal to be used for testing materials that frequently changes state from pressure to release and vice versa.

**Features:**

- Very simple operation
- Compact structure that is easy to maintain
- Excellent performance

**Technical specification**

Capacity of Test Specimen 30mm x 20mm

Frequency
Material
Voltage
Size of SS Plates
Number of Specimen Tested at a Time
Minimum Load



Astra Scientific

## Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265