

**Product Name :**  
Bench Thickness Gauge

**Product Code :**  
RUBBER SUK-AS03LAB006



**Description :**

Bench Thickness Gauge

**Technical Specification :**

**The Bench**

Thickness Gauge for rubber consists of a dial micrometer fitted on a rigid metallic frame. The micrometer has a flat circular indenter, which exerts the specified pressure through a dead weight on the rubber sample or test slab whose thickness is being measured.

**The**

rubber sample or slab is kept on an anvil, which is in form of a flat circular platform. The axis of the indenter and the anvil coincide and is the same as the axis of movement of the indenter. The indenter can be raised and lowered gently onto the rubber sample or slab with the help of a lifting lever.

The body of the thickness gauge is finished in shore blue metallic painting / black stove enamel painting. The anvil and the indenter are made from magnetic stainless steel. The loading weight is made from brass.

#### Technical specification

Anvil Diameter
Measuring Range
Material
Stainless Steel
Display
Resolution of Dial Micrometer
0.01mm

Register Load

Center Diameter



## Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265