

Product Name : Flex Testing Machine	Product Code : RUBBER SUK-AS03LAB009
---	--



Description :

Flex Testing Machine

Technical Specification :

The ability of shoes or chappals to withstand repeated bending experience during walking without developing a crack within their estimated life span is a very important factor for determining the quality of its sole. Also, the ability of the sole to withstand propagation of any slit formed due to a sharp edge or pointed object during walking is equally important in assessing its quality.

Technical specification

Bend Radius

Age

Frequency

Test Specimen Diameter

Width

Flexing Cycles Per Minutes

Bending Specimen Angle



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265