

Product Name :
Flex Testing Machine

Product Code :
RUBBER SUK-AS03LAB009



Description :

Flex Testing Machine

Technical Specification :

The ability of shoes or chappals to withstand repeated bending experience during walking without developing a crack within their estimated life span is a very important factor for determining the quality of its sole. Also, the ability of the sole to withstand propagation of any slit formed due to a sharp edge or pointed object during walking is equally important in assessing its quality.

Technical specification

--

Bend Radius

Wtage

50 Hzency

Test Specimen Diameter

315 Width

Fixing Cycles Per Minutes

Bending Specimen Angle



Astra Scientific

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265