

**Product Name :**  
Din Abrasion Tester

**Product Code :**  
RUBBER SUK-AS03LAB010



**Description :**

Din Abrasion Tester

**Technical Specification :**

- IS 3400 (Part 3) 1987 Methods of Test for Vulcanized Rubbers.
- Part 3 Abrasion Resistance using a Rotating Cylindrical Drum Device.

Technical Data

Diameter of test specimen
Height of test specimen

50 mm Diameter of rotating drum
470 gm Weight of rotating drum
800 RPM Speed of rotation of cylinder



Astra Scientific

## Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265