

Product Name :
Divider Flat Leg

Product Code :
LABUK-AP12LAB0001



Description :

Divider Flat Leg

Technical Specification :

- Made with forged steel legs.
- Hardened points for precise marking of the work piece.
- Quick adjustment.

- 10 inch flat legs and a 10-1/4 inch jaw capacity.
- Hardened and tempered bow springs and fulcrum spools to assure uniform tension without side play while using the divider.
- A flat leg divider is a useful measurement transfer tool for uneven surfaces.
- Also great for a variety of uses such as tracing contours and scribing circles into metal



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265