

Product Name :
Conductivity Meter

Product Code :
LABUK-AP12LAB0002



Description :

Conductivity Meter

Technical Specification :

- Conductivity Range:

1. 0 to 199.9 μ S,

2. 200 to 1999 μ S,

3. 2.00 to 19.99mS

• TDS/ Salinity Range :

1. 0 to 99.9 ppm (mg/ L)

2. 100 to 999ppm (mg/ L,

3. 1.00 to 9. 99ppt (g/L)

• pH Range: 0.00 to 14.00pH

• Temperature Range: 32° to 149°F (0° to 65°C)

• Max. Resolution: 0.1uS, 0.1ppm (mg/ L), 0.01pH, 0.1 °F/0c

• Accuracy: $\pm 2\%FS$, ± 0.01 pH, $\pm 1.8^\circ F/1^\circ C$

• Power: four SR44W button batteries

• Size: 1.4x7.3x1.6" (36x186x41mm)



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265