


**Product Name :**  
Magnetic Stirrer with Heating and Ceramic Heating Plate

**Product Code :**  
LABUK-AP12LAB0003



**Description :**

• Magnetic Stirrer with Heating and Ceramic Heating Plate

**Technical Specification :**

- Programmable Temp. Range: 5° C - 550°C.
- Type of Temperature Control: Digital.
- Over Temperature Protection: Yes

- Stirring Capacity: 15L of H2O
- Programmable Speed Range: 50- 1500 rpm
- Type Of Motor Control: Electronic
- Noise Level: < < 80dBa
- Temperature Range: +5 ... + 40 °C
- Storage Temperature Range: -10 ... +60 °C
- Maximum Humidity: 80%
- Maximum Altitude: 4000m
- Power Supply: 115V 60Hz
- Dimensions Width x Height x Depth : 203 X 94 X344mm
- Weight: 3.5 kg

- Overall Power: 800W
  
- Dimensions of the Heating Plate: 180 x180 mm.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265