

**Product Name :**  
Analytical Electronic Balance

**Product Code :**  
LABUK-AP12LAB0004



**Description :**

Analytical Electronic Balance

**Technical Specification :**

- Capacity: 220 g.

- Readability: (d) 0.1 mg
- Repeatability: 0.1 mg
- Linearity: +/- 0.3 mg
- Stabilization time: 3.0 Sec
- Ability to adjust the response time
- Pan size: 80 mm diameter
- Calibration self-calibration according to a temperature changes or pre- set conditions, would have Internal calibration
- Dimensions width x depth x height mm: 300x200x310
- Warm Up time: No , Plug and Play
- Wind shield parts would be easy to assemble and disassemble for cleaning purposes.
- Back-lit function for easy viewing.

- RS 232C port would be there to connect to printer or computer
- Power supply 220-230V 50 -60 Hz
- Measurement kit, Battery pack
- Warranty for 2 Years.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265