

**Product Name :**  
Vacuumed Chamber

**Product Code :**  
LABUK-AP12LAB0009



**Description :**

Vacuumed Chamber

**Technical Specification :**

Rough vacuum to 25 inches. (760 - 1 torr.)

Capacity: 75 liters

CFM/L/M displacement: 4

230VAC 50 Hz.

Warranty for 3 years.

Manufacture in India.

Operation manual provided.



## Astra Scientific

[www.astrascientific.com](http://www.astrascientific.com), **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265