

Product Name :
Electric Heating Mantle

Product Code :
LABUK-AP12LAB0010



Description :

Electric Heating Mantle

Technical Specification :

Features:

- Quick Heat-up time
- For round-bottom flasks of 1000 ml

- With rotary switch
- With an integrated controller.
- Maximum operation temperature: 0... +450 °C
- Voltage Supply 220 V.
- Power: 350W
- Operation time continuous Dimensions (mm): ø280 x 220
- Weight: 4 kg.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265