

Product Name :
Electric Oven

Product Code :
LABUK-AP12LAB0012



Description :

Electric Oven

Technical Specification :

Capacity: 75liter.

Shelf Positions: 3

No. Of Shelves: 2

Max Load Height (mm): 248 mm

Max Temperature (oC): 0 – 250oc

Fluctuation (oC): :±0.75

Power (Volts): 230 VAC

Power (Watts): 1250 W

Chamber: Aluminum Plated Steel

Power: 230 VAC /50Hz

Warranty for 3 years.

Manufacture in India.

Operation manual provided.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265