


Product Name :
Double Beam UV-VIS Spectrophotometer

Product Code :
LABUK-AP12LAB0018



Description :

Double Beam UV-VIS Spectrophotometer

Technical Specification :

Features:

- Type: Digital.
- Optical system: Double beam
- Wave length range : 190 to 1100 nm
- Wave length repeatability: 0.2 nm
- Spectrum bandwidth: 2.0 nm
- Photometric range : -0.3 to 3.0 A, 0 to 200%T, 0 to 9999.9 conc.
- Stray light: 0.1 %T
- Cell Holder: 4 position 10 mm path length.
- Photometric Accuracy : $\pm 0.5 T$
- Display mode : computer screen
- Wave length accuracy: $\pm 0.8 \text{ nm}$

- Wave length Resolution: 0.1nm
- Transmittance: 0 to 199T
- Absorbance: 0 to 1.000A
- Accuracy: $\pm 1\%T$
- Resolution: 0.1%T, 0.001A
- Mounting : Bench top
- Light source : Tungsten – Halogen, Deuterium lamp
- Display: LED display
- Power: Ac 230V , 50 Hz.

Accessories :

- Dust cover,

- 08(10 mm path length) glass curvets.
- 04 (10 mm path length)
- Quarts cuvetts.
- Warranty for 3 years.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265