

Product Name :
Digital Thermometer

Product Code :
LABUK-AP12LAB0021



Description :

Digital Thermometer

Technical Specification :

Features:

- Range: -50.0 to 150.0oC.
- Resolution: 0.1° C
- Accuracy: $\pm 0.5^{\circ}\text{C}$ (-20 to 90°C)
- Probe: Stainless Steel
- Switch: ON/OFF
- Calibration Heck.
- Environment: 0to 50°C; max RH 95%
- Battery Type: 1 x 1.5V (included)
- Battery Life 3000 hours of continuous
- Dimensions 66 x 50 x 25 mm
- Warranty: 2 years.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265