

Product Name :
Digital Clamp (Clip on) Meter

Product Code :
LABUK-AP12LAB0023



Description :

Digital Clamp (Clip on) Meter

Technical Specification :

Max flow rate: 1.5 l/s

3 - 3/4 Digit LCD with maximum reading of 4000

Overload Indication

Automatic Polarity

Low Battery Indication

ACV: 0- 750V

ACA:0-1200A

DC V: 0 - 1000 V

DC A: 0 - 1200 V

Min /Max.. Readings

Resistance: 0.1? to 40M ?

Continuity Buzzer

Diode testing

capacitance: 1pF to 30 pF

Frequency: 0.0Hz to 500 KHz Bar graph Auto ranging Auto Power off

Warranty for 3 years

Provided with operation manual



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265