

Product Name :
Digital Pressure Display (Precision Digital Pressure Gauge)

Product Code :
LABUK-AP12LAB0024



Description :

Digital Pressure Display (Precision Digital Pressure Gauge)

Technical Specification :

Measuring ranges up to 0-10 Bar vacuum and absolute pressure ranges

Accuracy: 0.025%

Intrinsically safe version

Fully temperature compensated accuracy

Panel mount gauges

Data logging option

Warranty for 3 years

Provided with operation manual



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265