


**Product Name :**  
Fluid Friction Apparatus

**Product Code :**  
LABUK-AP12LAB0025



**Description :**

Fluid Friction Apparatus

**Technical Specification :**

AISI 304 stainless steel structure

Total number of tapping points : 36

AISI 304 stainless steel pipes

Internal diameter of test pipes: For Different Sizes

Internal diameter rough pipe:

suitable diameter

Length of test pipe: 1000 mm

AISI 304 stainless steel 90o elbow

AISI 304 stainless steel 45o elbow

AISI 304 stainless steel 45o "Y"

AISI 304 stainless steel "T"

AISI 304 stainless steel sudden enlargement

AISI 304 stainless steel sudden contraction

Gate valve

Globe valve

"Y" strainer

Transparent venturimeter

AISI 304 stainless steel orifice plate

AISI 304 stainless steel Pitot tube

Water differential pressure gauge,

Range 0-1000 mm

Mercury differential pressure gauge

Range 0-1000 mm

Required services compatible

Hydraulic bench



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265