

Product Name :
Heavy Duty Oven

Product Code :
LABUK-AP12LAB0030



Description :

Heavy Duty Oven

Technical Specification :

Features:

- Capacity: 80Liter

- Shelves: 3

- Electric supply: 230 VAC, 50Hz, single
- Temperature Range : up to 200oC
- Control Accuracy:+ 2
- With Timer
- Resistance and durable finish
- Wire mesh trays
- Warranty: Minimum 1 year.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265