

|  |   |
|--|---|
| <b>Product Name :</b><br>Scrub Station Two Bay | <b>Product Code :</b><br>HOSPITALSUK-AS280002 |
|--|---|



|  |
|--|
| <b>Description :</b>   |
| Scrub Station Two Bay  |
| <b>Technical Specification :</b>                                     |
| With integrated sensor   |
| Digital timer with user defined time interval 1, 3, 5 and 10 minutes |
| Foot/knee operated manual mechanism                                  |

Metal body made of high grade stainless steel 304 grade

Provided with front access panel for ease of maintenance and repair



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265