

Product Name :
Scrub Station Two Bay

Product Code :
HOSPITALSUK-AS280002



Description :

Scrub Station Two Bay

Technical Specification :

With integrated sensor

Digital timer with user defined time interval 1, 3, 5 and 10 minutes

Foot/knee operated manual mechanism

Metal body made of high grade stainless steel 304 grade

Provided with front access panel for ease of maintenance and repair



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265