

Product Name :
X-Ray View Box LED with Dimmer and Sensor

Product Code :
HOSPITALSUK-AS330003



Description :

X-Ray View Box LED with Dimmer and Sensor

Technical Specification :

Used to view MRI, CT Scan and X Ray images

Energy efficient

Cost effective

Light Weight

Portable



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265