

Product Name :
Astra Automated Nucleic Acid Purification

Product Code :
ASTRALAB-SUPPLYA1015



Description :

Astra Automated Nucleic Acid Purification

Technical Specification :

There are two types of Auto-Pure series, Auto-Pure16A and Auto-Pure 32A. They are desktop devices. Purify DNA and RNA from different starting materials (such as blood, cultured cells or bacteria, tissues, cell-free body fluids and plant samples) by using magnetic separation technology.

The device can use a magnetic bar to transfer particles in each purification stage of binding, mixing, washing

and elution, thus providing a solution with the shortest hands-on time.

DNA and RNA purification has high quality and high molecular weight.

Auto-pure16A can accept up to 16 samples at the same time, while Auto-Pure 32A can accept up to 20 samples at the same time. The maximum sample size is 1ml.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265